

# Comparison sheet - Photizo Home Care units



| Product Name                   | Vetcare   | Pain Relief   | Sport   | Skin Care   | Biush   |
|--------------------------------|---|---|---|---|---|
| Diode Technology               | Light Emitting Diode (LED)  |
| CE Certification               | CE: For Animal use  | Medical CE: Chronic Pain  | Medical CE: Non-Chronic Pain  | Medical CE: Wound Healing   | Medical CE: Beauty Application  |
| Total Power Output             | 1300 mW   |
| Power Density (Peak / Average) | 340m W/cm2<br>194mW/cm2   |
| Wavelengths                    | 633nm (visible red)<br>850nm (invisible NIR)<br>both delivered simultaneously |
| Pulse / Continuous             | 30% pulsed - 70% continuous   | 30% pulsed - 70% continuous   | 30% pulsed - 70% continuous   | 30% pulsed - 70% continuous   | 30% pulsed - 70% continuous   |
| Timed Dose Duration            | 31 seconds  | 41 seconds  | 31 seconds  | 16 seconds  | 16 seconds  |
| Dose delivered                 | 6 joules/cm2  | 8 joules/cm2  | 6 joules/cm2  | 3 joules/cm2  | 3 joules/cm2  |
| Perform up to 75.000 hours     | ✓   | ✓   | ✓   | ✓   | ✓   |
| Charging Time                  | 3-4 hrs   |
| Mains Power                    | 100 – 240VAC  |
| Weight                         | 280g  | 280g  | 280g  | 280g  | 280g  |

For more information on the different units and how to treat pain, injury or wounds please don't hesitate to contact us on: [post@photizo.info](mailto:post@photizo.info)

